



ARL-TN-0948 • MAY 2019



Analysis of Pyrolysis Products of a Commercial Phenolic Resin by Gas Chromatography/Mass Spectrometry

by Rose Pesce-Rodriguez and James Wolbert

Approved for public release; distribution is unlimited.

NOTICES

Disclaimers

The findings in this report are not to be construed as an official Department of the Army position unless so designated by other authorized documents.

Citation of manufacturer's or trade names does not constitute an official endorsement or approval of the use thereof.

Destroy this report when it is no longer needed. Do not return it to the originator.



Analysis of Pyrolysis Products of a Commercial Phenolic Resin by Gas Chromatography/Mass Spectrometry

by Rose Pesce-Rodriguez and James Wolbert

Weapons and Material Research Directorate, CCDC Army Research Laboratory

REPORT DOCUMENTATION PAGE				Form Approved OMB No. 0704-0188	
<p>Public reporting burden for this collection of information is estimated to average 1 hour per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing the burden, to Department of Defense, Washington Headquarters Services, Directorate for Information Operations and Reports (0704-0188), 1215 Jefferson Davis Highway, Suite 1204, Arlington, VA 22202-4302. Respondents should be aware that notwithstanding any other provision of law, no person shall be subject to any penalty for failing to comply with a collection of information if it does not display a currently valid OMB control number.</p> <p>PLEASE DO NOT RETURN YOUR FORM TO THE ABOVE ADDRESS.</p>					
1. REPORT DATE (DD-MM-YYYY) May 2019		2. REPORT TYPE Technical Note		3. DATES COVERED (From - To) 3 April 2019–10 April 2019	
4. TITLE AND SUBTITLE Analysis of Pyrolysis Products of a Commercial Phenolic Resin by Gas Chromatography/Mass Spectrometry				5a. CONTRACT NUMBER	
				5b. GRANT NUMBER	
				5c. PROGRAM ELEMENT NUMBER	
6. AUTHOR(S) Rose Pesce-Rodriguez and James Wolbert				5d. PROJECT NUMBER	
				5e. TASK NUMBER	
				5f. WORK UNIT NUMBER	
7. PERFORMING ORGANIZATION NAME(S) AND ADDRESS(ES) US Army Combat Capabilities Development Command Army Research Laboratory ATTN: FCDD-RLW-LB Aberdeen Proving Ground, MD 21005.				8. PERFORMING ORGANIZATION REPORT NUMBER ARL-TN-0948	
9. SPONSORING/MONITORING AGENCY NAME(S) AND ADDRESS(ES)				10. SPONSOR/MONITOR'S ACRONYM(S)	
				11. SPONSOR/MONITOR'S REPORT NUMBER(S)	
12. DISTRIBUTION/AVAILABILITY STATEMENT Approved for public release; distribution is unlimited.					
13. SUPPLEMENTARY NOTES					
14. ABSTRACT A commercial phenolic resin was analyzed by pyrolysis-gas chromatography/mass spectrometry. Estimates of pyrolysis products (aromatic-phenolic, aromatic-nonphenolic, and permanent gas products) as a function of pyrolysis temperature are provided.					
15. SUBJECT TERMS phenolic resin, aromatic, methane, pyrolysis, gas chromatography/mass spectrometry, GC/MS, Cellobond J2027L					
16. SECURITY CLASSIFICATION OF:			17. LIMITATION OF ABSTRACT UU	18. NUMBER OF PAGES 53	19a. NAME OF RESPONSIBLE PERSON Rose Pesce-Rodriguez
a. REPORT Unclassified	b. ABSTRACT Unclassified	c. THIS PAGE Unclassified			19b. TELEPHONE NUMBER (Include area code) (410) 306-1877

Contents

List of Figures	iv
1. Introduction/Background	1
2. Experimental	1
3. Results and Conclusions	1
4. References	6
Appendix A. Mass Spectra for 250 °C Pyrolysis	7
Appendix B. Mass Spectra for 550 °C Pyrolysis	9
Appendix C. Mass Spectra for 650 °C Pyrolysis	13
Appendix D. Mass Spectra for 750 °C Pyrolysis	22
Appendix E. Mass Spectra for 850 °C Pyrolysis	34
Appendix F. Mass Spectra for 950 °C Pyrolysis	39
Appendix G. Mass Spectra for 1050 °C Pyrolysis	41
List of Symbols, Abbreviations, and Acronyms	43
Distribution List	44

List of Figures

Fig. 1	Total ion chromatogram for phenolic resin when pyrolyzed at 250 °C. Mass spectra given in Appendix A.	2
Fig. 2	Total ion chromatogram for phenolic resin when pyrolyzed at 550 °C	2
Fig. 3	Total ion chromatogram for phenolic resin when pyrolyzed at 450 °C	2
Fig. 4	Total ion chromatogram for phenolic resin when pyrolyzed at 550 °C. Assignments for peak times given in the insert. Mass spectra given in Appendix B.	3
Fig. 5	Total ion chromatogram for phenolic resin when pyrolyzed at 650 °C. Assignments for peak times given in the insert. Phenolic shown in black font, non-phenolic aromatics shown in gray font. Mass spectra given in Appendix C.	3
Fig. 6	Total ion chromatogram for phenolic resin when pyrolyzed at 750 °C. Assignments for peak times given in the insert. Phenolic shown in black font, non-phenolic aromatics shown in gray font. Mass spectra given in Appendix.	3
Fig. 7	Total ion chromatogram for phenolic resin when pyrolyzed at 850 °C. Assignments for peak times given in the insert. Phenolic shown in black font, non-phenolic aromatics shown in gray font. Mass spectra given in Appendix E.	4
Fig. 8	Total ion chromatogram for phenolic resin when pyrolyzed at 950 °C. Mass spectra given in Appendix F.	4
Fig. 9	Total ion chromatogram for phenolic resin when pyrolyzed at 1050 °C. Mass spectra given in Appendix G.	4
Fig. 10	Relative level of pyrolysis products platted against pyrolysis temperature. Permanent gas products are CO ₂ , CO, CH ₄ , and H ₂ O.	5
Fig. A-1	Mass spectrum of permanent gas peak from 250 °C pyrolysis of phenolic resin	8
Fig. B-1	Mass spectrum and library search match (with structure) of 3.0-min peak from 550 °C pyrolysis of phenolic resin	10
Fig. B-2	Mass spectrum and library search match (with structure) of 3.4-min peak from 550 °C pyrolysis of phenolic resin	10
Fig. B-3	Mass spectrum and library search match (with structure) of 3.5-min peak from 550 °C pyrolysis of phenolic resin	11
Fig. B-4	Mass spectrum and library search match (with structure) of 3.7-min peak from 550 °C pyrolysis of phenolic resin	11
Fig. B-5	Mass spectrum and library search match (with structure) of 3.9-min peak from 550 °C pyrolysis of phenolic resin	12

Fig. B-6	Mass spectrum and library search match (with structure) of 4.2-min peak from 550 °C pyrolysis of phenolic resin	12
Fig. C-1	Mass spectrum of permanent gas peak from 650 °C pyrolysis of phenolic resin	14
Fig. C-2	Mass spectrum and library search match (with structure) of 1.8-min peak from 650 °C pyrolysis of phenolic resin	14
Fig. C-3	Mass spectrum and library search match (with structure) of 2.0-min peak from 650 °C pyrolysis of phenolic resin	15
Fig. C-4	Mass spectrum and library search match (with structure) of 2.4-min peak from 650 °C pyrolysis of phenolic resin	15
Fig. C-5	Mass spectrum and library search match (with structure) of 3.4-min peak from 650 °C pyrolysis of phenolic resin	16
Fig. C-6	Mass spectrum and library search match (with structure) of 3.6-min peak from 650 °C pyrolysis of phenolic resin	16
Fig. C-7	Mass spectrum and library search match (with structure) of 3.8-min peak from 650 °C pyrolysis of phenolic resin	17
Fig. C-8	Mass spectrum and library search match (with structure) of 4.1-min peak from 650 °C pyrolysis of phenolic resin	17
Fig. C-9	Mass spectrum and library search match (with structure) of 5.7-min peak from 650 °C pyrolysis of phenolic resin	18
Fig. C-10	Mass spectrum and library search match (with structure) of 6.0-min peak from 650 °C pyrolysis of phenolic resin.	18
Fig. C-11	Mass spectrum and library search match (with structure) of 6.2-min peak from 650 °C pyrolysis of phenolic resin	19
Fig. C-12	Mass spectrum and library search match (with structure) of 6.3-min peak from 650 °C pyrolysis of phenolic resin	19
Fig. C-13	Mass spectrum and library search match (with structure) of 6.5-min peak from 650 °C pyrolysis of phenolic resin	20
Fig. C-14	Mass spectrum and library search match (with structure) of 6.7-min peak from 650 °C pyrolysis of phenolic resin	20
Fig. C-15	Mass spectrum and library search match (with structure) of 7.2-min peak from 650 °C pyrolysis of phenolic resin	21
Fig. D-1	Mass spectrum of permanent gas peak from 750 °C pyrolysis of phenolic resin	23
Fig. D-2	Mass spectrum and library search match (with structure) of 1.8-min peak from 750 °C pyrolysis of phenolic resin	23
Fig. D-3	Mass spectrum and library search match (with structure) of 2.0-min peak from 750 °C pyrolysis of phenolic resin	24
Fig. D-4	Mass spectrum and library search match (with structure) of 2.3-min peak from 750 °C pyrolysis of phenolic resin	24

Fig. D-5	Mass spectrum and library search match (with structure) of 2.9-min peak from 750 °C pyrolysis of phenolic resin	25
Fig. D-6	Mass spectrum and library search match (with structure) of 3.2-min peak from 750 °C pyrolysis of phenolic resin	25
Fig. D-7	Mass spectrum and library search match (with structure) of 3.4-min peak from 750 °C pyrolysis of phenolic resin	26
Fig. D-8	Mass spectrum and library search match (with structure) of 3.6-min peak from 750 °C pyrolysis of phenolic resin	26
Fig. D-9	Mass spectrum and library search match (with structure) of 3.9-min peak from 750 °C pyrolysis of phenolic resin	27
Fig. D-10	Mass spectrum and library search match (with structure) of 5.2-min peak from 750 °C pyrolysis of phenolic resin	27
Fig. D-11	Mass spectrum and library search match (with structure) of 5.5-min peak from 750 °C pyrolysis of phenolic resin	28
Fig. D-12	Mass spectrum and library search match (with structure) of 5.6-min peak from 750 °C pyrolysis of phenolic resin	28
Fig. D-13	Mass spectrum and library search match (with structure) of 5.9-min peak from 750 °C pyrolysis of phenolic resin.	29
Fig. D-14	Mass spectrum and library search match (with structure) of 6.0-min peak from 750 °C pyrolysis of phenolic resin	29
Fig. D-15	Mass spectrum and library search match (with structure) of 6.1-min peak from 750 °C pyrolysis of phenolic resin	30
Fig. D-16	Mass spectrum and library search match (with structure) of 6.2-min peak from 750 °C pyrolysis of phenolic resin	30
Fig. D-17	Mass spectrum and library search match (with structure) of 6.4-min peak from 750 °C pyrolysis of phenolic resin	31
Fig. D-18	Mass spectrum and library search match (with structure) of 6.5-min peak from 750 °C pyrolysis of phenolic resin	31
Fig. D-19	Mass spectrum and library search match (with structure) of 6.7-min peak from 750 °C pyrolysis of phenolic resin	32
Fig. D-20	Mass spectrum and library search match (with structure) of 7.1-min peak from 750 °C pyrolysis of phenolic resin	32
Fig. D-21	Mass spectrum and library search match (with structure) of 7.5-min peak from 750 °C pyrolysis of phenolic resin	33
Fig. E-1	Mass spectrum of permanent gas peak from 850 °C pyrolysis of phenolic resin	35
Fig. E-2	Mass spectrum and library search match (with structure) of 0.9-min peak from 850 °C pyrolysis of phenolic resin	35
Fig. E-3	Mass spectrum and library search match (with structure) of 1.0-min peak from 850 °C pyrolysis of phenolic resin	36

Fig. E-4	Mass spectrum and library search match (with structure) of 1.4-min peak from 850 °C pyrolysis of phenolic resin	36
Fig. E-5	Mass spectrum and library search match (with structure) of 2.2-min peak from 850 °C pyrolysis of phenolic resin	37
Fig. E-6	Mass spectrum and library search match (with structure) of 2.7-min peak from 850 °C pyrolysis of phenolic resin	37
Fig. E-7	Mass spectrum and library search match (with structure) of 6.6-min peak from 850 °C pyrolysis of phenolic resin	38
Fig. E-8	Mass spectrum and library search match (with structure) of 7.0-min peak from 850 °C pyrolysis of phenolic resin	38
Fig. F-1	Mass spectrum of permanent gas peak from 950 °C pyrolysis of phenolic resin	40
Fig. G-1	Mass spectrum of permanent gas peak from 1050 °C pyrolysis of phenolic resin	42

1. Introduction/Background

A sample of a commercial phenolic resin (i.e., Cellobond J2027L) was provided along with a request for an analysis of the pyrolysis products as a function of temperature. Specifically, there was an interest in the levels of permanent gas products (carbon monoxide [CO], carbon dioxide [CO₂], methane [CH₄], and water [H₂O]) as well as larger products (aromatic-phenolic and aromatic-nonphenolic) generated by pyrolysis of the resin from 200 to 1100 °C.

2. Experimental

Desorption and pyrolysis products were analyzed by means of a gas chromatography/mass spectrometry (GC/MS) instrument with a desorption interface. Desorption was achieved via a CDS Analytical Model 2000 Pyroprobe (coil type) connected through a heated interface chamber to the splitless injector of an Agilent (Santa Clara, California) GC/MS system (Model 6890N GC and Model 5973N MSD). The GC column used was a HP-5 capillary column (0.25 mm × 30 m, 0.25-μm film). The injector temperature was 200 °C; the Pyroprobe interface was set to a temperature of 250 °C. The GC oven temperature program was as follows: 100 °C isothermal for 1 min, 100–250 °C at 40 °C/min, and 250 °C isothermal for 1 min. The Pyroprobe was programmed to give a 20-s desorption pulses at a heating rate of 1,000 °C/s. The pulse temperature is based on calibration provided by the vendor and was not measured for this study. Samples (1–2 mg), which were prepared by curing Cellobond J2027L (Hexion Select) at 121 °C for 3 h in air, were held within the coil of the Pyroprobe by first placing them in a quartz tube containing a small plug of glass wool and then inserting the entire tube into the coil. Selected ion chromatograms were obtained via Hewlett Packard ChemStation software by extracting specified masses from the total ion chromatogram.

3. Results and Conclusions

Total ion chromatograms for pyrolysis of the phenolic resin are given in Figs. 1–9. Mass spectra, including best chemical identifications based on library matches, associated with each chromatogram are given in the Appendix indicated in the figure caption.




Fig. 1 Total ion chromatogram for phenolic resin when pyrolyzed at 250 °C. Mass spectra given in Appendix A.




Fig. 2 Total ion chromatogram for phenolic resin when pyrolyzed at 550 °C




Fig. 3 Total ion chromatogram for phenolic resin when pyrolyzed at 450 °C




Fig. 4 Total ion chromatogram for phenolic resin when pyrolyzed at 550 °C. Assignments for peak times given in the insert. Mass spectra given in Appendix B.




Fig. 5 Total ion chromatogram for phenolic resin when pyrolyzed at 650 °C. Assignments for peak times given in the insert. Phenolic shown in black font, non-phenolic aromatics shown in gray font. Mass spectra given in Appendix C.




Fig. 6 Total ion chromatogram for phenolic resin when pyrolyzed at 750 °C. Assignments for peak times given in the insert. Phenolic shown in black font, non-phenolic aromatics shown in gray font. Mass spectra given in Appendix.




Fig. 7 Total ion chromatogram for phenolic resin when pyrolyzed at 850 °C. Assignments for peak times given in the insert. Phenolic shown in black font, non-phenolic aromatics shown in gray font. Mass spectra given in Appendix E.




Fig. 8 Total ion chromatogram for phenolic resin when pyrolyzed at 950 °C. Mass spectra given in Appendix F.




Fig. 9 Total ion chromatogram for phenolic resin when pyrolyzed at 1050 °C. Mass spectra given in Appendix G.

Figure 10 gives a rough estimate of the levels of pyrolysis products generated at each temperature. The results are in general agreement with those in Katzman et al. (1995) and Bouajila et al. (2003), except for a somewhat higher temperature at which phenolic products are generated (i.e., approximately 700 °C in the current work versus approximately 500 °C the Katzman and Bouajila studies). However, the higher temperature is consistent with phenolic pyrolysis results from other works, including Chang and Tackett (1990), Szymański et al. (2002), and Wong et al. (2015), suggesting that comparisons cannot be made without carefully comparing the resin type and its structure, which is not possible with proprietary commercial products. For the given system in this study (i.e., Cellobond J2027L), it is concluded that pyrolysis occurs as indicated in the experimental results presented in this report. It is difficult to interpret the results further without knowing the specific composition of the commercial resin. Such an analysis is beyond the scope of this very limited investigation. Nevertheless, the pyrolysis GC/MS results herein provide valuable information on the chemical stability of Cellobond J2027L and on the types of gaseous decomposition that contribute to void formations at specific temperatures from 200 to 1100 °C.




Fig. 10 Relative level of pyrolysis products platted against pyrolysis temperature. Permanent gas products are CO₂, CO, CH₄, and H₂O.

4. References

- Bouajila J, Raffin G, Alamercury S, Waton H, Sanglar C, Grenier-Loustalot MF. Phenolic resins (IV). Thermal degradation of crosslinked resins in controlled atmospheres. *Polymers & Polymer Composites*. 2003;11(5). <http://www.polymerjournals.com/pdfdownload/894770.pdf>.
- Chang C, Tackett JR. Characterization of phenolic resins with thermogravimetry-mass spectrometry. Proceedings of the 19th Annual North American Thermal Analysis Society; 1990 Jan. OSTI Report No.: MLM-3668(OP); CONF-9009223-4; ON: DE91002710. <https://www.osti.gov/servlets/purl/6333297>.
- Katzman HA, Mallon JJ, Barry WT. Polyarylacetylene-matrix composites for solid rocket motor components. El Segundo (CA): Space and Missile Systems Center Air Force Materiel Command; 1995 Sep. Report No.: SMC-TR-95-48. DOI: <http://citeseerx.ist.psu.edu/viewdoc/download?doi=10.1.1.909.6606&rep=rep1&type=pdf>.
- Szymański GS, Karpiński Z, Biniak S, Świątkowski A. The effect of the gradual thermal decomposition of surface oxygen species on the chemical and catalytic properties of oxidized activated carbon. *Carbon*. 2002;40:2627–2639. DOI: [https://doi.org/10.1016/S0008-6223\(02\)00188-4](https://doi.org/10.1016/S0008-6223(02)00188-4).
- Wong H-W, Peck J, Bonomi RE, Assif J, Panerai F, Reinisch G, Lachaud J, Mansour NN. Quantitative determination of species production from phenol-formaldehyde resin pyrolysis. *Polymer Degradation Stability*. 2015;112:122–131. DOI: <https://doi.org/10.1016/j.polymdegradstab.2014.12.020>.

Appendix A. Mass Spectra for 250 °C Pyrolysis




Fig. A-1 Mass spectrum of permanent gas peak from 250 °C pyrolysis of phenolic resin

Appendix B. Mass Spectra for 550 °C Pyrolysis




Fig. B-1 Mass spectrum and library search match (with structure) of 3.0-min peak from 550 °C pyrolysis of phenolic resin




Fig. B-2 Mass spectrum and library search match (with structure) of 3.4-min peak from 550 °C pyrolysis of phenolic resin




Fig. B-3 Mass spectrum and library search match (with structure) of 3.5-min peak from 550 °C pyrolysis of phenolic resin




Fig. B-4 Mass spectrum and library search match (with structure) of 3.7-min peak from 550 °C pyrolysis of phenolic resin




Fig. B-5 Mass spectrum and library search match (with structure) of 3.9-min peak from 550 °C pyrolysis of phenolic resin




Fig. B-6 Mass spectrum and library search match (with structure) of 4.2-min peak from 550 °C pyrolysis of phenolic resin

Appendix C. Mass Spectra for 650 °C Pyrolysis




Fig. C-1 Mass spectrum of permanent gas peak from 650 °C pyrolysis of phenolic resin




Fig. C-2 Mass spectrum and library search match (with structure) of 1.8-min peak from 650 °C pyrolysis of phenolic resin




Fig. C-3 Mass spectrum and library search match (with structure) of 2.0-min peak from 650 °C pyrolysis of phenolic resin




Fig. C-4 Mass spectrum and library search match (with structure) of 2.4-min peak from 650 °C pyrolysis of phenolic resin




Fig. C-5 Mass spectrum and library search match (with structure) of 3.4-min peak from 650 °C pyrolysis of phenolic resin




Fig. C-6 Mass spectrum and library search match (with structure) of 3.6-min peak from 650 °C pyrolysis of phenolic resin




Fig. C-7 Mass spectrum and library search match (with structure) of 3.8-min peak from 650 °C pyrolysis of phenolic resin




Fig. C-8 Mass spectrum and library search match (with structure) of 4.1-min peak from 650 °C pyrolysis of phenolic resin




Fig. C-9 Mass spectrum and library search match (with structure) of 5.7-min peak from 650 °C pyrolysis of phenolic resin




Fig. C-10 Mass spectrum and library search match (with structure) of 6.0-min peak from 650 °C pyrolysis of phenolic resin.




Fig. C-11 Mass spectrum and library search match (with structure) of 6.2-min peak from 650 °C pyrolysis of phenolic resin




Fig. C-12 Mass spectrum and library search match (with structure) of 6.3-min peak from 650 °C pyrolysis of phenolic resin




Fig. C-13 Mass spectrum and library search match (with structure) of 6.5-min peak from 650 °C pyrolysis of phenolic resin




Fig. C-14 Mass spectrum and library search match (with structure) of 6.7-min peak from 650°C pyrolysis of phenolic resin




Fig. C-15 Mass spectrum and library search match (with structure) of 7.2-min peak from 650 °C pyrolysis of phenolic resin

Appendix D. Mass Spectra for 750 °C Pyrolysis




Fig. D-1 Mass spectrum of permanent gas peak from 750 °C pyrolysis of phenolic resin




Fig. D-2 Mass spectrum and library search match (with structure) of 1.8-min peak from 750 °C pyrolysis of phenolic resin




Fig. D-3 Mass spectrum and library search match (with structure) of 2.0-min peak from 750 °C pyrolysis of phenolic resin




Fig. D-4 Mass spectrum and library search match (with structure) of 2.3-min peak from 750 °C pyrolysis of phenolic resin




Fig. D-5 Mass spectrum and library search match (with structure) of 2.9-min peak from 750 °C pyrolysis of phenolic resin




Fig. D-6 Mass spectrum and library search match (with structure) of 3.2-min peak from 750 °C pyrolysis of phenolic resin




Fig. D-7 Mass spectrum and library search match (with structure) of 3.4-min peak from 750 °C pyrolysis of phenolic resin




Fig. D-8 Mass spectrum and library search match (with structure) of 3.6-min peak from 750 °C pyrolysis of phenolic resin




Fig. D-9 Mass spectrum and library search match (with structure) of 3.9-min peak from 750 °C pyrolysis of phenolic resin




Fig. D-10 Mass spectrum and library search match (with structure) of 5.2-min peak from 750 °C pyrolysis of phenolic resin




Fig. D-11 Mass spectrum and library search match (with structure) of 5.5-min peak from 750 °C pyrolysis of phenolic resin




Fig. D-12 Mass spectrum and library search match (with structure) of 5.6-min peak from 750 °C pyrolysis of phenolic resin




Fig. D-13 Mass spectrum and library search match (with structure) of 5.9-min peak from 750 °C pyrolysis of phenolic resin.




Fig. D-14 Mass spectrum and library search match (with structure) of 6.0-min peak from 750 °C pyrolysis of phenolic resin




Fig. D-15 Mass spectrum and library search match (with structure) of 6.1-min peak from 750 °C pyrolysis of phenolic resin




Fig. D-16 Mass spectrum and library search match (with structure) of 6.2-min peak from 750 °C pyrolysis of phenolic resin




Fig. D-17 Mass spectrum and library search match (with structure) of 6.4-min peak from 750 °C pyrolysis of phenolic resin




Fig. D-18 Mass spectrum and library search match (with structure) of 6.5-min peak from 750 °C pyrolysis of phenolic resin




Fig. D-19 Mass spectrum and library search match (with structure) of 6.7-min peak from 750 °C pyrolysis of phenolic resin




Fig. D-20 Mass spectrum and library search match (with structure) of 7.1-min peak from 750 °C pyrolysis of phenolic resin




Fig. D-21 Mass spectrum and library search match (with structure) of 7.5-min peak from 750 °C pyrolysis of phenolic resin

Appendix E. Mass Spectra for 850 °C Pyrolysis




Fig. E-1 Mass spectrum of permanent gas peak from 850 °C pyrolysis of phenolic resin




Fig. E-2 Mass spectrum and library search match (with structure) of 0.9-min peak from 850 °C pyrolysis of phenolic resin




Fig. E-3 Mass spectrum and library search match (with structure) of 1.0-min peak from 850 °C pyrolysis of phenolic resin




Fig. E-4 Mass spectrum and library search match (with structure) of 1.4-min peak from 850 °C pyrolysis of phenolic resin




Fig. E-5 Mass spectrum and library search match (with structure) of 2.2-min peak from 850 °C pyrolysis of phenolic resin




Fig. E-6 Mass spectrum and library search match (with structure) of 2.7-min peak from 850 °C pyrolysis of phenolic resin




Fig. E-7 Mass spectrum and library search match (with structure) of 6.6-min peak from 850 °C pyrolysis of phenolic resin




Fig. E-8 Mass spectrum and library search match (with structure) of 7.0-min peak from 850 °C pyrolysis of phenolic resin

Appendix F. Mass Spectra for 950 °C Pyrolysis




Fig. F-1 Mass spectrum of permanent gas peak from 950 °C pyrolysis of phenolic resin

Appendix G. Mass Spectra for 1050 °C Pyrolysis




Fig. G-1 Mass spectrum of permanent gas peak from 1050 °C pyrolysis of phenolic resin

List of Symbols, Abbreviations, and Acronyms

CCDC ARL	US Army Combat Capabilities Development Command Army Research Laboratory
CO	carbon monoxide
CO ₂	carbon dioxide
CH ₄	methane
GC/MS	gas chromatography/mass spectrometry
H ₂ O	water

1 DEFENSE TECHNICAL
(PDF) INFORMATION CTR
DTIC OCA

2 CCDC ARL
(PDF) IMAL HRA
RECORDS MGMT
FCDD RLD CL
TECH LIB

1 GOVT PRINTG OFC
(PDF) A MALHOTRA

2 CCDC ARL
(PDF) FCDD RLW LB
R PESCE-RODRIGUEZ
FCDD RLW MA
J WOLBERT