

Approved for public release; distribution unlimited.

Citation of trade names in this report does not constitute an official indorsement or approval of the use of such items.

Destroy this report when no longer needed. Do not return it to the originator.

12:13 Approved for public release; distriction unlimited. Citation of tilde names in this feased does not constitute an official indoreement of approval of the use of such items. 13.5 Destroy this report when no honger (needed. bo not return it to the originator.

Approved for public release; distribution unlimited.

.0.

1. 10 to 1.

(25 St.

AD 159124

TECHNICAL REPORT

73-26-FL

GLYCEROL INTERFERENCE IN THE DETERMINATION OF CRUDE FIBER. FAT, AND MOISTURE IN INTERMEDIATE MOISTURE MEAT PRODUCTS

By

Lloyd Cox



Project reference: 1J662713A034 Reproduced by NATIONAL TECHNICAL INFORMATION SERVICE 's (reportment of Commerce Springfreid VA 22151

Series: FL-159

August 1972

Food Laboratory US Army Natick Laboratories Natick, Mas.chusetts 01760

DOCUMEN	T CONTROL DATA -	R&D	lessified)
RIGINATING ACTIVITY (Corporate auth or)	uv ang annorenion müst	24. REPORT SECURITY CLASSIF	ICATION
US ARMY NATICK LABORATORIES		UNCLASSIFIED	
Natick, MA 01760		ZD. GROUP	
EPORT TITLE	<u></u>		
Glycerol Interference in the Det	ermination of Cr	ude Fiber,	
Fat, and Moisture in Intermediat	e Moisture Meat	Products	
DESCRIPTIVE NOTES (Type of report and inclusive dates))		
NO THORED (Part name, manale initial, tast name)			
Lloyd Cox			
EPORT DATE	Te TOTAL NO	OF PAGES 76 NO. OF REE	5
August 1972	-21	- 20 4	•
CONTRACT OR GRANT NO.	P. ORIGINAT	DR'S REPORT NUMBER(S)	*****
PROJECT NO. 1.THA271 24024	זיז	159	
	1.1.		
	Bb. OTHER P	NO(S) (Any other numbers that i	nay be assign
DISTRIBUTION STATEMENT			
SUPPLEMENTARY NOTES	US Ar Natic	ng Military activity my Natick Laboratories k, MA 01/760	
The AOAC crude fiber method f application to intermediate moistu fiber gave extremely high results, range. The AOAC Mojonnier (mixed food was modified and compared to Soxhlet method (petroleum ether) f higher fat values than the other t carbohydratcs. The Karl Fischer m and the AOAC toluene distillation Fischer results were comparable to moisture methods were compared for were comparable to the second day	for animal feed w are meat products the modified me ether) method for the Mojonnier (p for meat; the Moj two methods but t noisture method w method for moist to the AOAC toluer results on diff results.	as modified and evaluat ; where the AOAC method thod gave values in the r fat in baked or expan etroleum ether only) an c.nier (mixed cther) me here was no evidence of as compared to two oven ure in animal feed; the e distillation results. erent days; the first d	ed for for cruc expected ded pet d the A02 thod gave extracted mathods Karl The for ay result

• •

* *

*

all and the second

منهر از ال

単にい

•

, i ,

ł. ,

Unclassified Security Classification

Post WT Rott WT Rott WT Evaluation 8 9 7 1 1 1 Chemical Methods 9 7 4 7 1 1 Food Froducts 9 9 8 7 1 1 1 Applications 9 9 8 7 1 1 1 Food Froducts 9 9 8 9 1 1 1 Analysis Samples 9 8 1 1 1 1 Determination 0 1 1 1 1 1 1 Clycerol 1 1 1 1 1 1 1
Evaluation 8 9 7 Chemical Tets 9 9 7 Applications 9 9 9 Food Froducts 9 9 9 Analysis Samples 9 9 Moisture 0 4 9 Determination 0 4 4 Floer 9 4 4 Ploer 9 6 1 Pates 6 1 1 Glycerol 1 6 1

FOREWORD

One of the responsibilities of the Food Chemistry Division's Analytical Group is development of new chemical methods of analysis, or modification of existing ones for new and experimental food products.

1. ---- W. . .

When the high glycerol content of intermediate moisture meat products caused interference problems in the determinations of crude fiber, fat, and moisture, this study was conducted to evaluate the applicability of related AOAC* chemical methods. The results showed that selection and modification of existing methods of chemical analysis of food eliminated the interference problem.

This work was performed under Project 1.1662713A034, Military Food Service and Subsistence Technology.

Grateful acknowledgement is accorded personnel of these Laboratories: PFC John R. Troy and Mr. Otto J. Stark for their assistance in the statistical analysis of the data obtained, and Mr. Albert S. Henick, Dr. Maxwell C. Brockmann, and Dr. J. Walter Giffee for their review of the manuscript.

*Official Methods of Analysis (1970) 11th Ed., Association of Official Analytical Chemists, Washington, D.C., 1970.

j-b

TABLE OF CONTENTS

B.F ...

	Fage NO.
Foreword	i
List of Tables	iii
Abstract	iv
Introduction	1
Experimental	1
Method for Crude Fiber	2
Procedure	2
Results and Discussion	2
Method for Fat (Modified Mojonnier)	3
Procedure	3
Results and Discussion	3
Method for Moisture (Karl Fischer)	4
Procedure	4
Results and Discussion	5
Conclusions	5
References	7

LIST OF TABLES

Page No.

Table 1.	Comparative Crude Fiber (CF) Results of Ham and Raisin Sauce	8
Table 2.	Comparative Fat Results (%) of Pork with Sweet and Sour Sauce	9
Table 3.	Comparative Moisture Results (%) of Pork with Sweet and Sour Sauce	10
Table 4.	Comparative Moisture Results (%) of Barbecue Beef	11
Table 5.	Comparative Moisture Results (%) of Ham and Raisin Sauce	12
Table 6.	Comparative Moisture Results (%) of Hong Kong Pork	13
Table 7.	Comparison of Moisture Results (%) on Different Days	14

ABSTRACT

The AOAC crude fiber method for animal feed was modified and evaluated for application to intermediate moisture meat products; where the AOAC method for crude fiber gave extremely high results, the modified method gave values in the expected range. The AOAC Mojonnier (mixed ether) method for fat in baked or expanded pet food was modified and compared to the Mojonnier (petroleum ether only) and the AOAC Soxhlet method (petroleum ether) for meat; the Mojonnier (mixed ether) method gave higher fat values than the other two methods but there was no evidence of extracted carbohydrates. The Karl Fischer moisture method was compared to two oven methods and the AOAC toluene distillation method for moisture in animal feed; the Karl Fischer results were comparable to the AOAC toluene distillation results. The four moisture methods were compared for results on different days; the first day results were comparable to the second day results.

INTRODUCTION

The continuous development of new food products for the combat soldier often requires modification of present food chemistry methods as well as the development of new methods to overcome interfering components and give reliable and meaningful test data.

The infusion of glycerol (at the level of about 30% of the product) into meat products to increase stability and shelf-life produces interference problems in the determinations of crude fiber, fat, and moisture. In the analysis of crude fiber, glycerol was removed from the sample before the crude fiber was determined. In the Mojonnier fat determination, glycerol and other ethyl ether-soluble components were removed by washing the mixed ether extracts with water. In the determination of moisture, the Karl Fischer method was compared to two oven methods and the toluene distillation method for its applicability to the intermediate moisture products.

This study was conducted to evaluate chemical methods for application to intermediate moisture meat products.

EXPERIMENTAL

The Association of Official Analytical Chemists (AOAC), Official Methods of Analysis, 11th Edition, Chapters: Meat and Meat Products and Animal Feeds (1) was followed, with modifications, in the development of these methods.

The intermediate moisture samples were furnished by the Animal Products Division, Food Laboratory, US Army Natick Laboratories.

The samples were blended in a Waring* blender (Model CB-6 with small jar adapter) and stored in jars with air-tight screw, caps.

Three intermediate moisture samples of ham and raisin sauce were analyzed for crude fiber by the AOAC method and the results compared to those obtained with a modification of the method.

Two intermediate moisture samples of pork with sweet and sour sauce were analyzed for fat by the modified AOAC Mojonnier method (mixed ether extracts washed in a second Mojonnier extraction tube with water) and compared to the modified Mojonnier method (petroleum ether only) and the AOAC Soxhlet method (petroleum ether) for meat.

*Mention of company or trade names in this report does not imply endorsement over others not named.

A variety of intermediate moisture samples were analyzed for moisture by the Karl Fischer method (2, 3) and compared to the AOAC Animal Feed moisture method modified (vacuum oven, 16 hours at 70°C), the AOAC Meat and Meat Products moisture method (mechanical convection oven, 16-18 hours at 100-102°C), and the AOAC Animal Feed toluene distillation moisture method.

The four moisture methods were subjected to comparison of results on different days.

Method for Crude Fiber

When meat product samples containing large amounts of glycerol (about 30%) gave unusually high crude fiber values, a method was developed to remove glycerol from the sample before the crude fiber was determined. The modified method was compared to the AOAC method (7.053–7.057).

Apparatus and Reagents

- (a) AOAC 7.054 and 7.055.
- (b) Glassware Centrifuge tubes, 250 ml; beakers, 600 ml; stirring rods, 7 inches long; mortar and pestle.
- (c) Centrifuge With head for 250 ml centrifuge tubes.
- (d) Steam bath,
- (e) Diethyl ether Anhydrous, ACS grade.
- (f) Ethyl alcohol 95%, USP grade.

Procedure

Weigh 5-10g sample into 250 ml centrifuge tube. Extract (mix with stirring rod) with 75 ml ether, centrifuge, and decant to remove fat. Repeat the procedure. (Some glycerol is also removed.) Warm residue on steam bath with stirring to expel some of the ether. Add 50 ml of hot distilled water and heat on steam bath 10.15 minutes; add 50 ml alcohol, stir, centrifuge, and decant. Make two water alcohol extractions of the residue to remove the glycerol. Transfer residue to 600 ml beaker using alcohol wash bottle. Evaporate alcohol on steam bath or decant alcohol from 600 ml beaker and dry residue in 100°C oven. (Residue should be dry and crunchy.) Break up residue with flat-end stirring rod or grind in mortar and pestle. Determine crude fiber as in AOAC 7.053-7.057.

Results and Discussion

The data in Table 1 compare the AOAC crude fiber method with the modified method. The AOAC method for crude fiber analysis gave results 18-29 times higher than what is considered normal for ham and raisin sauce whereas the modified method provided values in the expected range.

The high values obtained from the AOAC crude fiber method were probably due to aldehyde crosslinking (tanning) in the samples by acrolein (4,5) when samples containing glycerol were dried (100-125°C in mechanical convection oven 1.2 hours) prior to removal of fat.

Method for Fat (Modified Mojonnier)

The modifications made to the AOAC Mojonnier method were: (1) larger sample size digested in beaker (2) elimination of use of alcohol (3) washing mixed ether extracts with water in a second Mojonnier extraction tube (flask) and (4) pouring the water-washed mixed ether extracts directly into a tared 150 ml beaker. The Mojonnier method, where only petroleum ether was used did not require the water washing step.

Larger samples (10g) were used in the Soxhlet method for fat (AOAC 24.005), and petroleum ether was used instead of ethyl ether.

Procedure

Weigh 3-5g into 50 ml beaker and make paste with water (2.4 ml). Add 8-10 ml of conc. HCI (digest sample with 8N HCI if cereal content is high) and mix immediately with stirring rod. Cover beaker with watch glass and digest sample on steam bath for 30-45 minutes with occasional stirring. Transfer digested sample into Mojonnier extraction tube with distilled water and add sufficient water to reach the pouring off level. Add 25 ml ethyl ether (in two 12.5 ml increments) into 50 ml beaker with ether wash bottle to extract traces of fat residue and transfer into Mojonnier extraction tube. Stopper Mojonnier tube with cork and mix gently. Carefully release pressure and wash cork with a few ml of petroleum ether using wash bottle. Add 25 ml of petroleum ether (in two 12.5 ml increments) into 50 ml beaker and transfer into Mojonnier tube. Stopper Mojonnier tube, shake, and centrifuge (600 rpm) for 3-5 minutes. Transfer mixed ether extract (use small funnel) into second Mojonnier tube containing water, stopper, shake, centrifuge, and transfer water washed extract into tared 150 ml beaker. (Crystallizing dish, 80 x 40 mm. may be used instead of 150 ml beaker.) Evaporate extract in forced-draft hood or on steam bath. Make three additional extractions by pouring 25 ml of each ether directly into the first Mojonnier tube and following the same sequence of steps. (All transfer steps are made in a quantitative manner.) After the last extract evaporates, dry extracted fat in a mechanical convection oven at 100°C to constant weight (1.5-2 hrs.), cool in desiccator, and weigh.

Results and Discussion

The data in Table 2 permit comparison of fat recovered by the three methods. The mixed ether Mojonnier method with water washing in a second Mojonnier extraction tube gave the highest recoveries of fat, but no charred material (evidence of extracted carbohydrate, etc.) was found in the extracted fat.

In a similar study by Kuhn (6), the fat recovered by the Mojonnier method was higher than fat recovered by the Soxhlet method (petroleum ether). The Kuhn study also showed that, when the residue from the Soxhlet method was subjected to the Mojonnier method and the percent fat recovered was added to the percent fat by the Soxhlet method, the two fat values (Mojonniei value compared to total value of Soxhlet plus additional Mojonnier of Soxhlet residue) agreed reasonably well.

The data in Table 2 were subjected to analysis of variance, and the F-values showed significant difference at 0.05 probability when the Mojonnier method using mixed ether was compared to the Mojonnier method using petroleum ether only or the Soxhlet method using petroleum ether.

Method for Moisture (Karl Fischer)

The Karl Fischer method for moisture (2,3) was applied to intermediate moisture meat samples.

Apparatus and Reagents

- (a) Karl Fischer Aquametry Apparatus Labindustries, Berkeley, CA 94710.
- (b) Oster blender With blending assembly for standard Mason jar.
- (c) Mason jar 1/2 pint.
- (d) Karl Fisher Reagent (KFR) Harleco 3786, Arthur H. Thomas Co., Philadelphia, PA 19105.
- (e) Methanol Anhydrous.
- (f) Standard Weight 1g H_2O , to nearest mg, into 100 ml volumetric flask and bring to volume with methanol. Use flask short enough to fit on analytical balance, and cap flask with aluminum foil while weighing to prevent loss of H_2O . Titrate 5 ml aliquot of standard to determine KFR equivalence (eq.).

Procedure

Weight 5g sample into 1/2 pint Mason jar; immediately add 95 ml of methanol, and cap jar. Attach Oster blending assembly to jar and blend sample three 10 second intervals Allow sample to settle in tightly capped jar. Blank for sample should be determined on 95 ml of methanol blended with dry blending assembly in the same manner as sample. Equilibrate Karl Fischer (KF) apparatus by adding sufficient methanol and a slight excess of KFR with stirring for about 15 minutes. Add enough 98% methanol (a few drops) to give a slight excess of H₂O, and titrate to exact KF end-point. Pipette 5 ml of standard, sample, or blank into KF cell, and titrate to KF end-point. Empty cell by aspiration. $%H_2O = (ml \text{ KFR} - \text{blank}) \times \text{KFR}$ eq. x 100/mg sample, where mg sample = 5g x 1000/100 ml x 5 ml = 250 mg.

Results and Discussion

10.21 2 5

The data in Tables 3-6 show that the Karl Fischer method compared favorably to the toluene distillation method, with the exception of the values in Table 6.

The values for the vacuum oven method were slightly higher than the values for the Karl Fischer method and the toluene distillation method. The higher values were due partly to evaporation of acetic acid from the sample. Entrainment (extremely fine drops carried away with the vapor) of glycerol from the samples was indicated by condensation on the vacuum oven door.

The values for the mechanical convection oven method were extremely high compared to the other three methods. The data suggest that entrainment of glycerol from the samples is greater at atmospheric pressure than at reduced pressure.

Newman (7) observed the problem of loss of glycerol on 24 hours .5 drying glycerol solutions at reduced pressure and recommended Karl Fischer moisture determination for complex mixtures containing glycerol. Lawrie (4) observed that when distilling dilute glycerol solutions, loss due to entrainment was much greater at atmospheric pressure than at reduced pressure.

The data in Tables 3-6 were subjected to analysis of variance. The F-values comparing the toluene distillation method to the Karl Fischer method indicated insignificant difference at probability 0.05 for each product except Hong Kong Pork (Table 6) where the mean value for the Karl Fischer method is about 0.8% higher than the mean value for the toluene distillation method. The F-values comparing the mechanical convection oven and the vacuum oven method to the Karl Fischer method indicated significant differences at probability 0.05 in each product tested.

In Table 7 the four moisture methods were compared for results on different days. The data were subjected to analysis of variance. The F values comparing the first day to the second day indicated insignificant differences in each method except the vacuum oven method where at 0.05 probability there was slight significant difference.

CONCLUSIONS

When fat and glycerol were removed from the intermediate moisture meat product as described in this report, low and meaningful test results were obtained for crude fiber.

The Mojonnier fat method (mixed ether) gave the highest test results, and the water washing eliminated interference of glycerol, carbuhydrates, etc. The method can probably be applied to semi-moist pet food products and other food products containing ethyl ether soluble components (glyce', glycerol, carbohydrates, lactic acid, sorbate, etc.) The Karl Fischer and toluene distillation methods were equally good for determination of moisture in intermediate moisture meat products. The vacuum oven and mechanical convection oven methods gave high results, especially the latter method. The vacuum oven method is not recommended, but it can be used if consideration is made for evaporation of acetic acid and entrainment of a small percent of glycerol from the sample.

REFERENCES

- (1) Official Methods of Analysis (1970) 11th Ed., AOAC, Washington, D.C.
- (2) Mitchell, J., Jr., and Smith, D. M. (1948) Aquametry, Interscience Publications, Inc., New York.
- (3) Labindustries Aquametry Apparatus Instructions, Labindustries, 1802 Second Street, Berkeley, CA.
- (4) Lawrie, J. W. (1928) Glycerol and the Glycols, The Chemical Catalog Co., Inc., New York, pp 56-61, 187.

•

×.* . *

:

- (5) Gustavson, K. H. (1956) The Chemistry of Tanning Processes, New York Academic Press, New York, pp. 244-245.
- (6) Kuhn, G. O., Collaborative Study of the Determination of Fat in Semi-moist Pet Foods, Quaker Oats Research Labs, Barrington, IL. (This paper was presented to the 85th Annual Meeting of the AOAC, 1971.)
- (7) Newman, A. A. (1968) Glycerol, C.R.C. Press, Cleveland, pp. 16, 20, 95-96.

Sample #	% CF, AOAC	% CF, MODIFIED AOAC
03259	6.72	0.36
	6.72	0.39
03260	5.47	0.28
	5.74	0.24
03261	5.71	0.20
	5.30	0.18

Table 1. Comparative Crude Fiber (CF) Results of Ham and Raisin Sauce

N

	Sa	ample #12217		S	ample #1 221 8	
	Mojonniër, Mixed Ether	Mojonnier, Pet Ether	Soxhlet, Pet Ether	Mojonnier, Mixed Ether	Mojonnier, Pet Ether	Soxhlet, Pet Ether
	5.15	4.82	4.83	5.59	5.31	5.29
	5.11	4.75	4.80	5.57	5.20	5.26
	4.95	4.86	4.82	5.43	5.27	5.29
	4.98	- 4.83	4.87	5.39	5.25	5.29
	5.01	4.89	4.83	5.38	5.32	5.26
Mean	5.04	4.83	4.83	5.47	5.27	5.27
Std. Dev.	0.09	0.05	0.02	0.10	0.05	0.01
F value		21.39	27.39		16.36	18.11

Table 2. Comparative Fat Results (%) of Pork with Sweet and Sour Sauce

T.

;

. . .

	Karl Fischer	Vacuum Oven	Mechanical Convection Oven	Toluene Distillation
	38.02	38.54	45.74	38.0
	37.90	38.78	35.12	38.0
	37.54	38.89	46.32	38.2
	37.78	39.53	46.39	37.8
	38.02	39.42	45.26	
	38.25	39.37	46.71	
Mean	37.91	39.10	45.92	38.0
Std. Dev.	0.24	0.40	0.65	0.16
F-value		36.93	799.0	0.34

Table 3. Comparative Moisture Results (%) ofPork with Sweet and Sour Sauce

	Karl Fischer	Vacuum Oven	Mechanical Convection Oven	Toluene Distillation
	42.03	43.45	48.60	41.0
	41.91	43.06	49.33	41.5
	42.26	43.61	47.36	42.0
		43.66	48.28	41.5
		43.85	48.81	
		43.85	47.89	
Mean	42.07	43.58	48.38	41.5
Std. Dev.	0.18	0.30	0.70	0.41
F-value		63.71	224.1	4.89

Table 4. Comparative Moisture Results (%) of Barbecue Beef

1.67

÷

2

. . . .

1 14 Name

١,

State R. S.

	Karl Fischer	Vacuum Oven	Mechanical Convection Oven	Toluene Distillation
	41.68	43.79	48.93	42.0
	41.79	43.37	48.96	41.8
	42.15	44.21	50.92	42.7
	41.91	44.08	51.86	42.5
	42.15	43.34	50.60	
		43.46	50.85	
		42.93	51.86	
		42.88	50.94	
Mean	41.94	43.51	50.62	42.3
Std. Dev.	0.21	0.49	1.13	0.42
F-value		45.02	279.81	2.16

Table 5. Comparative Moisture Results (%) ofHam and Raisin Sauce

142479月27

	Karl Fischer	Vacuum Oven	Mechanical Convection Oven	Toluene Distillation
	42.62	42.90	48.12	42.0
	42.50	42.89	47.70	41.5
	42.62	43.00	50.60	42.0
	42.62	43.14	51.02	41.5
	42.59	43.04	50.70	41.5
	43.17	43.15	50.01	42.0
	42.82	43.27	50.58	41.5
		43.54	49.55	42.0
				41.5
······································				42.0
Mean	42.59	43.12	49.79	41.8
Std. Dev.	0.06	0.21	1.25	0.29
F-value		22.18	126.72	32.52

Table 6. Comparative Moisture Results (%) of H ng Kong Purk

ALL HALL

16-1-14 24-13-20-201

and the second s

;

• .

Table 7. Comparison of Moisture Results (%) on Different Days

-

		Pork	with Sweet a	nd Sour Sauce			Hong Kon	ig Pork
	Karl	Fischer	Vacuun	n Oven	Mechan Convectio	ical n Oven	Tolue Distillati	ane ion
	1st day	2nd day	1st day	2nd day	1st day	2nd day	1st day	2nd day
	37.66	37.90	38.54	39.81	45.74	44.06	42.0	41.5
	38.02	37.78	38.78	39.33	45.12	44.46	41.5	42.0
	37.54	38.25	38.89	39.11	46.32	44.17	42.0	41.5
	37.42		38.53	38.99	46.39	46.08	A1.5	42.Û
	38.02		39.37	40.02	45.26	45.88		41.5
			39.42	39.69	\$6.71	45.31		42.0
Mean	37.73	37.98	38.92	39.49	45.92	44.99	41.8	41.8
Std. Dev.	0.28	0.24	0.39	0.41	0.65	0.88	0.29	0.27
F-values		1.59	9	.05	4	32	0.0	00

:

• :

14

.

* *****

•••