

# REPORT DOCUMENTATION PAGE

*Form Approved*  
OMB No. 0704-0188

Public reporting burden for this collection of information is estimated to average 1 hour per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden to Washington Headquarters Services, Directorate for Information Operations and Reports, 1215 Jefferson Davis Highway, Suite 1204, Arlington, VA 22202-4302, and to the Office of Management and Budget, Paperwork Reduction Project (0704-0188), Washington, DC 20503.

<b>1. AGENCY USE ONLY (Leave blank)</b>	<b>2. REPORT DATE</b> 15 January 1999	<b>3. REPORT TYPE AND DATES COVERED</b> Final 1 October 1997-30 September 1998
---	--	---

<b>4. TITLE AND SUBTITLE</b> Children's Hospital of Pittsburgh Histocompatibility Center Molecular Typing of Alleles at HLA Class I and Class II Loci on Blood Spotted Filter Paper	<b>5. FUNDING NUMBERS</b> N00014-97-1-1060
--	---

<b>6. AUTHOR(S)</b> Massimo Trucco, M.D.	
---	--

<b>7. PERFORMING ORGANIZATION NAMES(S) AND ADDRESS(ES)</b> Children's Hospital of Pittsburgh Rangos Research Center 3460 Fifth Avenue Pittsburgh, PA 15213-3205	<b>8. PERFORMING ORGANIZATION REPORT NUMBER</b> 335
---	--

<b>9. SPONSORING / MONITORING AGENCY NAMES(S) AND ADDRESS(ES)</b> Department of the Navy Office of Naval Research 800 North Quincy Street Arlington, VA 22217-5660	<b>10. SPONSORING / MONITORING AGENCY REPORT NUMBER</b>
--	---

19990120 105

<b>11. SUPPLEMENTARY NOTES</b>
--------------------------------

<b>a. DISTRIBUTION / AVAILABILITY STATEMENT</b> Approved for Public Release; distribution is Unlimited.	<b>12. DISTRIBUTION CODE</b>
--	------------------------------

**13. ABSTRACT (Maximum 200 words)**

The activities of the Children's Hospital of Pittsburgh Histocompatibility Center can be classified under three major categories: 1) the repository service; 2) the molecular typing service; and 3) the research activities aimed at improving the technical support on which the previous two services are working. This year in the Repository, where blood samples from volunteers as bone marrow donors are collected from more than 120 recruiting centers nationwide, roughly 150,000 samples were processed. Close to 450,000 vials were stored in three sets of freezers (one sample per donor each freezer), while nearly 130,000 were sent out to other Typing Laboratories. In our laboratories, approximately 10,000 typings were performed with a very high Q.C. standard, as monitored by National Marrow Donor Program. Our Center's research branch synthesized all the oligonucleotides to be used as PCR primers or labelled probes to perform molecular HLA typing. DNA sequencing was also performed when hybridization results were not optimal. Alternatives for storing blood samples (e.g., blood spots on filter paper) or for enhancing typing capabilities (e.g., microchip technology to expedite PCR procedures) were successfully implemented or improved to the point of justifying their future implementation.

<b>14. SUBJECT TERMS</b>	<b>15. NUMBER OF PAGES</b> 7
	<b>16. PRICE CODE</b>

<b>17. SECURITY CLASSIFICATION OF REPORT</b> Unclassified	<b>18. SECURITY CLASSIFICATION OF THIS PAGE</b> Unclassified	<b>19. SECURITY CLASSIFICATION OF ABSTRACT</b> Unclassified	<b>20. LIMITATION OF ABSTRACT</b> UL
--	---	--	---

CHILDREN'S HOSPITAL OF PITTSBURGH HISTOCOMPATIBILITY CENTER

AWARD NUMBER: N00014-97-1-1060

Total Samples Typed and Results Sent from 10/01/97 to 09/30/98

	Total	Priority 1	Priority 2	Priority 3	No Makes	Navy
<b>1997</b>						
October	1013	119	894	0	0	0
November	851	161	690	0	0	0
December	822	122	700	0	3	0
<b>1998</b>						
January	838	130	708	0	0	0
February	854	263	591	0	1	0
March	1029	348	681	0	2	0
April	820	174	646	0	0	0
May	976	116	860	0	0	0
June	865	300	565	0	0	0
July	805	49	756	0	0	0
August	818	159	659	0	0	0
September	1026	188	838	0	0	0

**CHILDREN'S HOSPITAL OF PITTSBURGH HISTOCOMPATIBILITY CENTER**

AWARD NUMBER: N00014-97-1-1060

Total Samples Received, Stored and Shipped from the Repository from 10/01/97 to 09/30/98

<b>SAMPLES</b>	<b>NUMBER (October 1997)</b>	<b>NUMBER (November 1997)</b>	<b>NUMBER (December 1997)</b>
Stored	11,593	13,743	12,523
Shipped (AB - Class I)	4,700	5,800	7,400
Shipped (DR - Class II)	9,394	6,614	8,463
Not stored	87	113	149
Destroyed	0	238	0
<b>SAMPLES</b>	<b>NUMBER (January '98)</b>	<b>NUMBER (February '98)</b>	<b>NUMBER (March '98)</b>
Stored	8,192	7,437	16,264
Shipped (AB - Class I)	3,400	4,300	5,100
Shipped (DR - Class II)	3,578	5,713	5,652
Not stored	47	28	124
Destroyed	1,736	1,941	0

<b>SAMPLES</b>	<b>NUMBER (April '98)</b>	<b>NUMBER (May '98)</b>	<b>NUMBER (June '98)</b>
Stored	14,769	14,028	11,628
Shipped (AB - Class I)	5,250	4,100	5,700
Shipped (DR - Class II)	4,862	4,580	4,140
Not stored	173	99	98
Destroyed	464	0	1,131
<b>SAMPLES</b>	<b>NUMBER (July '98)</b>	<b>NUMBER (August '98)</b>	<b>NUMBER (September '98)</b>
Stored	11,709	10,703	10,601
Shipped (AB - Class I)	5,100	4,300	3,175
Shipped (DR - Class II)	4,031	4,668	5,115
Not stored	101	143	59
Destroyed	1,594	588	2,039

- Stored samples (Multiply number stored x 3 because of triplicate aliquots)
- Shipped samples (Separated by AB and DR samples)
- Not Stored because of insufficient quantity; duplicate tubes; too clotted; broken; wrong anti-coagulant.
- Destroyed because no longer eligible, so that blood sample must be removed from NMDP list.

CHILDREN'S HOSPITAL OF PITTSBURGH HISTOCOMPATIBILITY CENTER

AWARD NUMBER: N00014-97-1-1060

Research Completed from 10/01/97 to 09/30/98

Month/ Year:	Oligonucleotides Synthesized:	ABI/PE Sequencing:	Visible Genetics Sequencing:	Research Projects:
Oct '97	31	163 reactions	0*	<p>Began exploring optimal PCR conditions for both HLA-A and HLA-B typing using Visible Genetics Microblaster Sequencer.</p> <p>Continued to expand cell reference lines for cell and DNA bank.</p> <p>Synthesized 1 Taqman Probe.</p>
Nov '97	26	42 reactions	0*	<p>New Rhodamine terminator chemistry validated on sequencer.</p> <p>Rhodamine matrix installed on sequencer.</p> <p>Installation of Visible Genetics Sequencer.</p> <p>Class 1 HLA-A and HLA-B PCR amplification optimized for sequencing conditions.</p>
Dec '97	16	104	14	<p>Performed control samples for Class 1 typing on Visible Genetics sequencer.</p> <p>Continued to optimize amplifications and cycle sequencing for HLA-A and HLA-B.</p>
Jan '98	24	94	5	<p>Began training personnel on Visible Genetics Sequencer for Class 1 typing.</p> <p>Performed a quality control study for Visible Genetics sequencer using 3 unknown typing samples and 3 known typing samples.</p>

Month/ Year:	Oligonucleotides Synthesized:	ABI/PE Sequencing:	Visible Genetics Sequencing:	Research Projects:
Feb '98	62	65	11	<p>Set up data base for Visible Genetics Class 1 typing.</p> <p>Continued validation of Class 1 typing via sequencing.</p> <p>Troubleshooting on Class 1 typing reactions using Visible Genetics sequencer.</p>
Mar '98	31	130	7	<p>Big Dye terminator chemistry validated and installed on sequencer.</p> <p>Optimized cycle sequencing for Big Dye Terminator Kit.</p> <p>Expanded cell line references.</p> <p>Continued troubleshooting on Class 1 typing reactions using Visible Genetics sequencer.</p>
Apr '98	26	88	2	<p>Continued to optimized cycle sequencing for Big Dye Terminator Kit.</p> <p>Evaluated and optimized PCR "clean-up" methods for new sequencing chemistry.</p> <p>Continued to expanded cell line references.</p>
May '98	19	117	3	<p>Installed and updated Visible Genetics sequencer : OpenStep 3.0 software.</p> <p>Synthesized 3 Taqman Probes.</p> <p>Began evaluating DNA sequencing systems using capillary electrophoresis on a computer chip.</p> <p>Continued to expanded cell line references.</p>

Month/Year:	Oligonucleotides Synthesized:	ABI/PE Sequencing:	Visible Genetics Sequencing:	Research Projects:
Jun '98	43	89	11	<p>Evaluated sequencing methods for chip technology.</p> <p>Developed chip design for DNA sequencing.</p> <p>Initiated orders for chip sequencing.</p> <p>Installed Iomega "Jazz" System tape back-up for the hard drive on the Visible Genetics Sequencer.</p> <p>Synthesized 1 Taqman Probe.</p> <p>Continued to expanded cell line references.</p>
Month/Year:	Oligonucleotides Synthesized:	ABI/PE Sequencing:	CE Sequencing**:	Research Projects:
July '98	41	124	12	<p>Began instrumental set-up for DNA sequencing via capillary electrophoresis (CE) - computer chip technology.</p> <p>Optimized PCR for CE sequencing</p>
Aug '98	28	86	18	<p>Explored polymers and began initial testing.</p> <p>Continued instrumental set-up for CE chip sequencing runs.</p> <p>Continued optimizing PCR conditions for CE runs.</p>
Sep '98	22	125	10	<p>Continued optimizing PCR conditions for CE runs.</p> <p>Continued to evaluate polymers for CE sequencing.</p> <p>Began literature search on new HPLC techniques for DNA purification.</p>

\* Visible Genetics sequencing instrument was not yet purchased at this time.

\*\* Instituted new project utilizing capillary electrophoresis (CE) for DNA sequencing in July 1998.