UNCLASSIFIED

AD NUMBER ADB175871 NEW LIMITATION CHANGE TO Approved for public release, distribution unlimited **FROM** Distribution authorized to U.S. Gov't. agencies only; Test and Evaluation; 28 JUL 1993. Other requests shall be referred to Army Medical Research and Development Command, ATTN: SGRD-RMI-S, Fort Detrick, MD 21702-5012. **AUTHORITY** AMRMC MCMR-RMI-S [70-1y] memo. dtd 28 Aug 1995



U.S. ARMY MEDICAL RESEARCH AND MATERIEL COMMAND FORT DETRICK, MARYLAND 21702-5012



EPLY TO ATTENTION OF MCMR - RMI - S

(70-1y)

MEMORANDUM FOR Admin Defense Technical Information Center, ATTN: DTIC-HDS/William Bush,

Cameron Station, Building 5, Alexandria, VA

22304-6145

SUBJECT: Request Change in Distribution Statement

The U.S. Army Medical Research and Materiel Command has reexamined the need for the limited distribution statement on technical report for Military Interdepartmental Purchase Request Number 92MM2556. Request the limited distribution statement for AD Number ADB175871 be changed to "Approved for mublic wedges; distribution unlimited." A copy of this report should be released to the National Technical Information Service.

2. Point of contact for this request is Ms. Judy Pawlus, DSN 343-7322.

COL, MS

Deputy Chief of Staff for Information Management

REPRODUCTION QUALITY NOTICE

This document is the best quality available. The copy furnished to DTIC contained pages that may have the following quality problems:

- Pages smaller or larger than normal.
- Pages with background color or light colored printing.
- Pages with small type or poor printing; and or
- Pages with continuous tone material or color photographs.

Due to various output media available these conditions may or may not cause poor legibility in the microfiche or hardcopy output you receive.

If this block is checked, the copy furnished to DTIC contained pages with color printing, that when reproduced in Black and White, may change detail of the original copy.

AD____

MIPR NO: 92MM2556

TITLE: PHYSIOLOGICAL STUDY OF 4 MODIFIED SELF-CONTAINED

VENTILATORS

SUBTITLE: Evaluation of the Impact 78G Self-Contained

Ventilator

PRINCIPAL INVESTIGATOR: Peter H. Abbrecht, M.D., Ph.D.

CONTRACTING ORGANIZATION: Uniformed Services University

of the Health Sciences
4301 Jones Bridge Road

Bethesda, Maryland 20814-4799

REPORT DATE: July 28, 1993

TYPE OF REPORT: Final Report

PREPARED FOR: U.S. Army Medical Research and

Development Command, Fort Detrick Frederick, Maryland 21702-5012

DISTRIBUTION STATEMENT: Distribution authorized to U.S. Government agencies only, Test and Evaluation, July 28, 1993. Other requests for this document shall be referred to Commander, U.S. Army Medical Research and Development Command, SGRD-RMI-S, Fort Detrick, Frederick, Maryland 21702-5012.

The findings in this report are not to be construed as an official Department of the Army position unless so designated by other authorized documents.

93-20342

93 8 31 062

(Z)

impaired pulmonary function. 14. SUBJECT TERMS 15. NUMBER OF PAGES Mechanical Ventilation; Ventilator; Resuscitation 16. PRICE COCE Device; RAD V SECURITY CLASSIFICATION 18. SECURITY CLASSIFICATION 19. SECURITY CLASSIFICATION 20. LIMITATION DE AH THAL OF REPORT OF THIS PAGE OF ABSTRACT Unclassified Unclassified Unclassified Limited

N.N 7540001 230 5500

Standard Form (144 5) Pinu bing by ANN (13) (18 8 198 (12) Project No: G176ER

Evaluation of the Impact 780
Self Contained Ventilator

Peter H. Abbrecht, M.D., Ph.D.

Uniformed Services University of the Health Sciences

4301 Jones Bridge Road Bethesda, MD 20814-4799

nusu . jul	UNENIATION PA		
	car continuating to payingge the set girls	·	Mark to the second seco
me. and maintaining the sate wedge and the	signing and reviewing the collection of it	nthreation was headed in doubted berides birectorate	garding this builden estimate or any illner assess in this coordinates of the control of the con
1. AGENCY USE ONLY ,Leave blank)			NO DATES COVERED
4. TITLE AND SUBTITLE	7-28-1993	Final	4-01-1992 12-31-199 Ts. FUNDING NUMBERS
S. 1112 Aug 300 1112			J. VONDING NORDERS
Evaluation of the Impac	t 780 Self-Contain	ed Ventilator	9216 42556
6. AUTHOR(S)			1
Peter H. Abbrecht, M.D	., Ph.D.		
7. PERFORMING ORGANIZATION NAME	(S) AND ADDRESS(ES)		8. PERFORMING ORGANIZATION
			REPORT NUMBER
Uniformed Services Univ 4301 Jones Bridge Rd Bethesda, MD 20814-4	•	th Sciences	
9. SPONSORING, MONITORING AGENC	NAME(S) AND ADDRESS(ES		10. SPONSORING MONITORING
U.S. Army Medical Mater Ft. Detrick Frederick, MD 21702-9	-	tivity	AGENCY REPORT NUMBER
12a. DISTRIBUTION / AVAILABILITY STAT	TEMENT		12b. DISTRIBUTION CODE
13. ABSTRACT (Maximum 200 words)			
before and after vibrate compressor and external impedances to air flow circuit, and the C-2 at including minute ventilitime, exhalation time, different operating moventilatory performance changed significantly corresponding change is leaked outward during inward during spontaneous breathing a spontaneous breathing a spontaneous inspiration.	tion testing. Unit compressed air por of the input and our filter, and determined the lation, tidal volume and duty cycle. We des. There was relie of the four differing response to vibration and the resistance of the was 26.5 cm H ₂ O/(L/mathrough the vent).	s were tested u wered ventilati utput ports, th rmined the vent e, mean deliver also determine atively little rent units. Th ation was insuf ll four units, filter connecto during powered four units to sec), a value t	on modes. We measured to spontaneous inspiration ilatory parameters and flow rate, insufflation distance in the difference in the conly variable that flation time, with a the exhalation ports are leaked unfiltered air eventilation with the inspiratory airflow during hat would probably make
impaired pulmonary fund		-	
14. SUBJECT TERMS			115. NUMBER OF PAGES

Mechanical Ventilation

16 PRICE CODE

17 SECURITY CLASSIFICATION OF THIS PAGE OF ABSTRACT

Unclassified Unclassified Unclassified CAP

NSN 7540-01 280-5500

Standard form 298 Her 2 ha

The Impact 780 Self-Contained Ventilator consists of an electrically powered compressor, non-adjustable control circuitry, internal battery, external electrical power converter, air filter, hoses, and a multi-valve assembly to join the patient to the ventilator. The internal compressor draws in filtered ambient air and delivers it to the patient. High pressure fittings and a hose are provided to use compressed gas (50 to 80 psi) as an alternate form of powered ventilation. Electrical power is required for either type of powered ventilation, and can be supplied by the external power converter and battery, the external power converter alone, or the battery alone. The patient can always inspire filtered air in addition to, or in place of, either type of powered ventilation. Patient spontaneous inspiration does not require electrical power.

Four ventilators were tested at USUHS to determine the values of the mission critical clinical parameters, both before and after vibration challenge (IAW MIL-STD-810 Vibration, Method 514.3, Procedure X, XI), with three different sources of electrical power. The ventilators were evaluated using the internal compressor and external compressed gas powered ventilation modes. All compressed-gas ventilation testings were performed with compressed air. We also determined the flow/pressure relationships and compliance that would be experienced by a patient breathing spontaneously through the ventilator.

1. General Observations

The exhalation ports on all four patient valves (connector fitting between the patient and the ventilator tubing) leaked outward during powered ventilation insufflation. All filter connectors leaked unfiltered air inward during spontaneous inspiration and during powered ventilation with the internal compressor. In the latter case, from 111 to 765 ml/min of unfiltered air leaked past the filter into the system and was delivered to the patient valve. This leakage around the filter would also expose the patient to unfiltered air during spontaneous inspiration. We were able to eliminate the filter/connector leaks with gaskets. One of the units also had an outward leak inside the ventilator on the positive pressure side of the internal compressor circuit (419 ml/min) which reduced the amount of ventilation delivered to the patient.

The controls and hose fittings on the ventilator were so close together that they were difficult to use. One fuze blew inside a ventilator and another in a power converter. Just to validate the status of the fuzes required the use of tools to open the units resulting in exposure of their delicate internal components. The cause(s) of the blown fuzes was never discovered. Once the fuzes were replaced, there were no further occurrences of this problem.

The external compressed-gas circuit has a low pressure electronic alarm and warning light and a high pressure whistle. The low pressure alarm was active at an average value of 45 psi or less before the vibration challenge, and 47 psi after. The high pressure alarm turned on at an average of 89 psi before the vibration challenge, and 95 psi after. Once activated the high pressure alarm then turned off when the pressure was reduced to an average of 76 psi or less before the vibration challenge, and 84 psi after.

When operating with battery power, battery life was shorter in the internal compressor mode than when using external compressed gas. In the internal compressor mode starting with a fully charged battery, the average time to the first, intermittent low

Dist Special

battery alarm buzzer and light was 5.6 hours before the vibration challenge, and 4.8 hours after. The average time to continuous low battery alarm was 6.6 hours before the vibration challenge, and 6.8 hours after. Compressor action stopped completely after an average 7.2 hours operation before the vibration challenge, and 7.9 hours after.

The patient valve has a pressure relief valve to prevent exposing the patient to excessive pressures. The relief valve opened at an average pressure of 31 cm H₂O before the vibration challenge, and 28 cm H₂O after. Escaping flow through the relief valve sounded an alarm at an average pressure of 36 cm H₂O before the vibration challenge, and 35 cm H₂O after. For pressures between the opening level and the alarm level, gas was vented without any audible indication.

During powered ventilation, the maximum pressure that the patient might be exposed to is determined by the amount of flow delivered to the patient valve and by the characteristics of the pressure relief valve. With the internal compressor powered by the external power converter, the unit average peak pressure inside the patient valve was 46 cm H₂O before the vibration challenge and 43 cm H₂O after. With the internal compressor powered by just the internal battery, the average maximum pressure inside the patient valve was 44 cm H₂O before the vibration challenge and 41 cm H₂O after. Using 50 psi external compressed gas, the average maximum patient valve internal pressure was 42 cm H₂O before the vibration challenge and 40 cm H₂O after, and was only slightly influenced by choice of electrical power source.

2. Ventilator Impedance

The amount of ventilation delivered to the patient during powered ventilation depends on the ventilator circuit resistance, the patient's pulmonary resistance, and the amount of power (electrical or gas pressure) driving the ventilator. The following are average resistances at 1 L/sec flow rate measured before and after the vibration challenge.

	Resistance, cm H ₂ O/(L/sec) ± S			
Vibration Challenge State	Before	After		
Patient Valve Components				
Powered Ventilation Input Port	231 ± 26	264 ± 37		
Spontaneous Inspiration Input Port	6.8 ± 0.3	6.3 ± 0.4		
Exhalation Output Port, without PEEP	7.5 ± 1.3	9.2 ± 2.5		
Exhalation Output Port, with PEEP	19.5 ± 1.9	20.2 ± 2.0		
Pressure Relief Valve	72.7 ± 2.6	79.3 ± 5.8		
Entire Spontaneous Inspiration Circuit	26.3 ± 1.6	26.7 ± 1.8		
Spontaneous Inspiration Hose	6.9 ± 0.5	7.4 ± 0.2		
C-2 Air Filter	2.4 ± 0.1	2.4 ± 0.1		

The average resistance (26.5 cm H₂O/(L/sec)) of the four units to inspiratory airflow during spontaneous breathing would probably make spontaneous inspiration through the ventilator impossible for a person with impaired pulmonary function.

To test the effects of increased patient impedance on ventilator performance, the systems were tested with seven different resistances at the patient interface port during powered ventilation using both the internal compressor and external compressed air. The following are the average tidal volumes during different powered ventilation modes for two extreme patient resistances.

	Tidal Volu	mes. L ± SD
Added Resistance, cm H ₂ O/(L/sec)	0	50
Internal Compressor, Internal Battery and External Power Converter	0.96 ± 0.06	0.82 ± 0.03
Internal Compressor, Internal Battery	0.82 ± 0.04	0.74 ± 0.03
External Compressed Air at 50 psi Internal Battery and External Power Converter	0.77 ± 0.02	0.72 ± 0.02

3. Ventilator Capability

Compressed-Gas Ventilation

The effect of different inlet pressures on tidal volumes during compressed-air ventilation was determined. The following are the average tidal volumes for supplied inlet air pressures spanning the inlet pressure alarm limits.

	Tidal Volumes, L					
Inlet Pressure, psi	50	60	70	80		
Tidal Volumes (liters)	0.76	0.90	1.03	1.18		

Spontaneous-Inspiration Compliance

The spontaneous inspiration circuit average static compliance was almost constant for a range of vacuums up to -25 cm H_2O . The average compliance was 0.41 ml/cm H_2O before the vibration challenge and 0.46 ml/cm H_2O after.

Powered Ventilation Performance

The following tables summarize the powered ventilation results for each of the two powered ventilation modes (internal compressor and external compressed air at 50 psi). Each ventilation mode was tested with three difference sources of electrical power (external power converter and battery, external power converter alone, and battery alone). For each ventilation mode and power source, 10 measurements were made on each of the four units. Thus, the mean values are the average of 40 measures (ten replications of each of the four units).

Powered Ventilation Mode: Internal Compressor Electrical Power Source: External Power Converter and Battery

Before Vibration Challenge

Ventilation Variables	(units)	Mean	Std Dev	Std Error
Minute Ventilation,	L/min	15.48	0.90	0.14
Tidal Volume,	L	0.978	0.054	0.008
Mean Delivered Flow Rate,	L/sec	0.571	0.034	0.005
Frequency,	1/min	15.83	0.32	0.05
Duty Cycle,	%	45.18	1.07	0.17
Insufflation Time,	sec	1.713	0.048	0.008
Expiration Time,	sec	2.078	0.064	0.010

Ventilation Variables	(units)	Mean	Std Dev	Std Error
Minute Ventilation,	L/min	15.07	0.86	0.14
Tidal Volume,	L	0.959	0.054	0.009
Mean Delivered Flow Rate,	L/sec	0.550	0.029	0.005
Frequency,	1/min	15.72	0.19	0.03
Duty Cycle,	9/0	45.72	0.68	0.11
Insufflation Time,	sec	1.746	0.031	0.005
Expiration Time,	sec	2.074	0 (39	Orman

Powered Ventilation Mode; Internal Compressor Electrical Power Source: External Power Converter

Before Vibration Challenge

Ventilation Variables	(units)	Mean	Std Dev	Std Error
Minute Ventilation,	L/min	15.15	0.92	0.15
Tidal Volume,	· L	0.960	0.056	0.009
Mean Delivered Flow Rate,	L/sec	0.561	0.037	0.006
Frequency,	1/min	15.79	0.34	0.05
Duty Cycle,	%	45.03	1.00	0.16
Insufflation Time,	sec	1.712	0.058	0.009
Expiration Time,	sec	2.090	0.052	0.008

Ventilation Variables	(units)	Mean	Std Dev	Std Error
Minute Ventilation,	L/min	14.71	0.74	0.12
Tidal Volume,	L	0.939	6.048	0.008
Mean Delivered Flow Rate,	L/sec	0.537	0.026	0.004
Frequency,	1/min	15.67	0.21	0.03
Duty Cycle,	%	45.74	0.72	0.12
Insufflation Time,	sec	1.752	0.033	0.005
Expiration Time,	sec	2.079	0.042	0.007

Powered Ventilation Mode: Internal Compressor Electrical Power Source: Battery

Before Vibration Challenge

Yentilation Variables	(units)	Mean	Std Dev	Std Error
Minute Ventilation,	L/min	12.89	0.73	0.12
Tidal Volume,	L	0.824	0.046	0.007
Mean Delivered Flow Rate,	L/sec	0.483	0.024	0.004
Frequency,	1/min	15.64	0.24	0.04
Duty Cycle,	%	44.55	1.13	0.18
Insufflation Time,	sec	1.710	0.045	0.007
Expiration Time,	sec	2.128	0.060	0.009

Ventilation Variables	(units)	Mean	Std Dev	Std Error
Minute Ventilation,	L/min	12.78	0.56	0.09
Tidal Volume,	L	0.817	0.036	0.006
Mean Delivered Flow Rate,	L/sec	0.473	0.019	0.003
Frequency,	1/min	15.64	0.23	0.04
Duty Cycle,	%	45.02	0.85	0.13
Insufflation Time,	sec	1.727	0.041	0.006
Expiration Time,	sec	2.110	0.046	0.007

Powered Ventilation Mode: External Compressed-Air (50 psi at inlet during flow) Electrical Power Source: External Power Converter and Battery

Before Vibration Challenge

Ventilation Variables	(units)	Mean	Std Dev	Std Error
Minute Ventilation,	L/min	12.21	0.49	80.0
Tidal Volume,	L	0.767	0.015	0.002
Mean Delivered Flow Rate,	L/sec	0.494	0.016	0.002
Frequency,	1/min	15.91	0.44	0.07
Duty Cycle,	%	41.22	0.79	0.13
Insufflation Time,	sec	1.555	0.034	0.005
Expiration Time,	sec	2.220	0.081	0.013

Ventilation Variables	(units)	Mean	Std Dev	Std Error
Minute Ventilation,	L/min	12.18	0.37	0.06
Tidal Volume,	L	0.775	0.021	0.003
Mean Delivered Flow Rate,	L/sec	0.495	0.018	0.003
Frequency,	1/min	15.72	0.20	0.03
Duty Cycle,	%	40.97	0.64	0.10
Insufflation Time,	sec	1.565	0.030	0.005
Expiration Time,	sec	2.255	0.039	0.006

Powered Ventilation Mode: External Compressed-Air (50 psi at inlet during flow) Electrical Power Source: External Power Converter

Eefore Vibration Challenge

Ventilation Variables	(units)	<u>Mean</u>	Std Dev	Std Error
Minute Ventilation,	L/min	12.12	0.39	0.06
Tidal Volume,	L	0.767	0.015	0.002
Mean Delivered Flow Rate,	L/sec	0.493	0.017	0.003
Frequency,	1/min	15.81	0.31	0.05
Duty Cycle,	%	40.95	0.77	0.12
Insufflation Time,	sec	1.555	0.042	0.007
Expiration Time,	sec	2.243	0.052	0.008

Ventilation Variables	(units)	Mean	Std Dev	Std Error
Minute Ventilation,	L/min	12.21	0.35	0.06
Tidal Volume,	L	0.778	0.020	0.003
Mean Delivered Flow Rate,	L/sec	0.497	0.017	0.003
Frequency,	1/min	15.68	0.20	0.03
Duty Cycle,	%	40.94	0.58	0.09
Insufflation Time,	sec	1.567	0.031	0.005
Expiration Time,	sec	2.261	0.035	0.006

Powered Ventilation Mode: External Compressed-Air (50 psi at inlet during flow) Electrical Power Source: Battery

Before Vibration Challenge

Ventilation Variables	(units)	Mean	Std Dev	Std Error	
Minute Ventilation,	L/min	12.04	0.37	0.06	
Tidal Volume,	L	0.767	0.0014	0.002	
Mean Delivered Flow Rate,	L/sec	0.491	0.014	0.002	
Frequency,	1/min	15.69	0.31	0.05	
Duty Cycle,	%	40.91	0.93	0.15	
Insufflation Time,	sec	1.566	0.039	0.006	
Expiration Time,	sec	2.262	0.064	0.010	

Ventilation Variables	(units)	Mean	Std Dev	Std Error
Minute Ventilation,	L/min	12.24	0.38	0.06
Tidal Volume,	L	0.781	0.022	5.003
Mean Delivered Flow Rate,	L/sec	0.498	0.017	0.003
Frequency,	1/min	15.67	0.22	0.04
Duty Cycle,	%	40.97	0.53	0.08
Insufflation Time,	sec	1.568	0.029	0.005
Expiration Time,	sec	2.261	0.038	0.006

The compressor motor used during internal compressor ventilation requires much more power from the battery than does the solenoid during external compressed-gas ventilation. The internal compressor produces a higher flow rate when supplied by the external power converter than when supplied by just the internal battery. Therefore, the internal compressor mode has larger variability for the minute ventilation, tidal volume, and mean delivered flow rate as a function of electrical power source than does the external compressed-gas mode.

The p values of the main effects from the analyses of variance are also listed. The analysis assessed the main effects and their interactions of three independent factors. These factors were: the units, the vibration challenge (VIB), and the sources of electrical power. There were 4 units of the same design. There were two VIB states: before and after. There were three different sources of electrical power: external electrical power converter and battery, external electrical power converter alone, and battery alone. The ten replicated measures were treated as nested terms. We chose p values less than or equal to 0.05 as statistically significant, and they are indicated in bold print.

The following is a listing of p values of the main effects from an analysis of variance for each dependent ventilation variable and each ventilation mode. Each analysis included ten replications for each repeated measure of four units, before and after the vibration challenge, using the three different electric power sources.

				p Values			
	Internal Compressor		External Compressed Air				
	Unit	VIB	Power		Unit	VIB	Power
Minute Ventilation,	0.0001	0.1133	0.0001		0.0001	0.6748	0.0027
Tidal Volume,	0.0001	0.2068	0.0001		0.0001	0.3516	0.1258
Mean Delivered Flow Rate,	0.0001	0.0539	0.0001		0.0001	0.5942	0.5605
Frequency,	0.0001	0.1331	0.0008		0.0001	0.1676	0.0017
Duty Cycle,	0.0001	0.0367	0.0045		0.0001	0.5780	0.1768
Insufflation Time,	0.0001	0.0100	0.0748		0.0001	0.0632	0.1183
Expiration Time,	0.0001	0.3571	0.0007		0.0001	0.2651	0.0045

While there were many statistically significant differences for the various ventilation variables, the differences may not be functionally significant especially because the variation in ventilation requirements of different sized patients may markedly exceed the small differences in ventilator performance observed in these tests.